



Low Cost 5mm NMR Sample Tubes Designed for Bench Top and High Field Rapid Scan Applications

Quality and Value you can rely on!®

Bench Top

Prices in US Dollars, valid through December 2021

Designed for **Nanalysis** Bench Top NMR instruments, manual load only.

4.95mm +/- 0.03mm OD x 0.43mm +/- 0.02mm Wall

Made of Type 1, Class B borosilicate glass.

Each pack of tubes comes with caps.

Request pricing for larger quantities.

NE-BT5-7-50	pack of 50 tubes, 7", with caps	\$54.95/pack
NE-BT5-7-100	pack of 100 tubes, 7", with caps	107.79/pack
NE-BT5-8-50	pack of 50 tubes, 8", with caps	58.12/pack
NE-BT5-8-100	pack of 100 tubes, 8", with caps	114.13/pack



High Field Rapid Scan

Designed for high field NMR instruments with rapid scan (high through-put) applications with/wo auto changer.

4.95mm +/- 0.02mm OD x 0.43mm +/- 0.02mm Wall

Made of Type 1, Class B borosilicate glass.

Each pack of tubes comes with caps.

Request pricing for larger quantities.

NE-RST5-7-50	pack of 50 tubes, 7", with caps	\$58.12 /pack
NE-RST5-7-100	pack of 100 tubes, 7", with caps	114.13/pack
NE-RST5-7-250	pack of 250 tubes, 7", with caps	274.76/pack
NE-RST5-8-50	pack of 50 tubes, 8", with caps	61.29/pack
NE-RST5-8-100	pack of 100 tubes, 8", with caps	120.47/pack
NE-RST5-8-250	pack of 250 tubes, 8", with caps	290.62/pack



Extra Caps

Standard 5mm tube cap, made of EVA (ethylene vinyl acetate)

NE-310-5-Orange-100	pack of 100 caps, Orange	\$9.57/pack
NE-310-5-Orange-1000	pack of 1000 caps, Orange	64.42/pack
NE-310-5-Purple-100	pack of 100 caps, Purple	9.57/pack
NE-301-5-Purple-1000	pack of 1000 caps, Purple	64.42/pack



Extra Labels

A clear label with a 10 x 19mm white writing area and a wrap-around feature to cover and protect the notation. **For spinning applications, position the label close to the spinner for the best results.**

NE-7408-BT-100	roll of 100 labels	\$33.82 / roll
----------------	--------------------	----------------

